Axis UL-8^{WSLA}



Powerchair

Freedom to move where you want when you want. Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.

- > Shoprider ultra-light powerchair
- > Easy as 1-2-3 to disassemble
- > 12AH batteries
- > Off-board 2A charger
- > Colour: red







UL-8^{WSLA} Specifications

Weight capacity	Dimensions LXWXH	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
113кg/250lbs	30"x 21.5"x 32"	47кg/104lbs	2 х .5нр	6крн/3.9мрн	2 х 12ан	f6" r9"	14км/9м	40см/15.5"	8°/14%
All specifications subject to change without notice.									

Dealer stamp



Discovering your world again-in style

www.eclipsemedical.com