Spyder XL P327 PATHMASTER

Powerchair

Eclipse Pathmaster Series – Freedom to move where you want when you want. Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and outdoor applications.

- Pathmaster Spyder XL 12 inch mid-wheel powerchair
- > Full suspension
- > Off-board 5A charger with 55AH batteries
- > Swing-away joystick
- > Colours: red or blue







P327 Specifications

Weight capacity	Dimensions LXWXH	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb ^o
205kg/450lbs	35"x 27"x 50"	106кg/233lbs	2 х .8нр	8крн/5мрн	2 х 55ан	F6" м12" r6"	32км/20м	53см/21"	8°/14%
All specifications subject to change without notice.									

Dealer stamp



Discovering your world again-in style

www.eclipsemedical.com