Navigator^L P424^L



Powerchair

Freedom to move where you want when you want.

Manoeuverability and affordability are the features offered by our line of power chairs. The full range of power chairs meet the individual needs of the consumer. Perfect for both indoor and

outdoor applications.



- > Full suspension
- > Off-board 5A charger with 50AH batteries
- > Swing-away joystick
- > Colours: red, blue, black or champagne with interchangeable panels included













$P424^{L}$ Specifications

Total wt. w/bat. Turning radius Safe climbo Weight capacity Dimensions LXWXH Motor **Speed** Driving dist. **Battery** 158kg/350lbs 41"x 25"x 46" 115kg/253lbs 2 х 1.34нр 9.6крн/6мрн 2 х 50ан F6" M14" R5" 34km/21m 57cm/22.5" 9°/15%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style